

AMERICAN BIOTECH LABS

SilverSol Technology®



On April 30 2009, American Biotech Labs held an informational webinar about its uniquely patented SilverSol Technology®. For your convenience, we have prepared this handout to summarize the information that was presented by Dr. Gordon Pedersen.

Webinar outline

- What is SilverSol Technology®
- Recommended uses of both liquid and gel
- Supporting University and doctor research
- Mechanisms of Action
- SilverSol® and the Flu Virus
- Case studies
- Before and after pictures

"...the Silver killed the residual bacteria that the antibiotic couldn't."

- *Current Science, 2006*

"In short, we currently do not have anything with such a wide spectrum of efficacy... in our inventory."

- *Paul K. Carlton, Jr., M.D., FACS*

SilverSol Technology®

- Approved in US Patent #7135195
- Destroys Bacteria, Viruses, and Fungi
- Does not harm good bacteria
- Can be taken every day for Prevention
- Can be used in every orifice of the body
- Helps reduce inflammation
- Permanently Suspended
- Resonates between 890-910 Terahertz (Frequency similar to UV light)
- Helps stimulate Stem Cells
- Used on body both internally and externally
- Helps regenerate wounds on the skin

SilverSol Technology® and Flu Virus

ABL has animal model viral studies on H5N1 as reported in a peer reviewed article published by JSHO which states:

"Especially of interest was the observation that 60% of the infected mice treated with this compound survived compared to the 30% in the placebo-treated controls"

Viral Studies:

In test tube results the American Biotech Labs silver product was able to kill or neutralize a billion bacteriophage viruses in 2.5 hours. In test tube tests against an Avian Influenza A H3-N2 and the Beijing Influenza A H1-N1, the American Biotech Labs 10 ppm

product was able to kill or neutralize 96-98% of the virus in two hours, with no measurable virus surviving in 12 hours. Against Avian Influenza A H5-N1 Vietnam Hybrid, the 10 ppm product reduced the viral levels below detectable levels in six hours.

Keith Moeller, a Managing Director states, "We realize that these study findings are just a preliminary step in showing the effectivity of the product as a daily supplement to improve health. But, given the fact that there is no other product that has been proven both safe for daily usage, and even moderately effective against H5-N1 Bird Flu, we think the study is very important."



SilverSol® can be taken daily as a prophylactic supplement



SilverSol[®] Liquid



SilverSol[®] Gel

SilverSol[®] Liquid

- Clear, odorless
- Virtually Tasteless
- No preservatives, additives, or artificial ingredients
- Used both orally, as well as topically

SilverSol[®] Gel

- Clear, odorless gel
- Contains no alcohol, or petroleum
- Improves wound healing
- Can be applied topically all over the body
- An extra-strength version just FDA approved

Recommended Doses	Recommended Doses
<ul style="list-style-type: none"> • Prevention <ul style="list-style-type: none"> ○ Swallow one teaspoon ○ 1-3 times a day • Post Infection Immune Support <ul style="list-style-type: none"> ○ Two teaspoons twice a day 	<ul style="list-style-type: none"> • Prevention <ul style="list-style-type: none"> ○ Apply to hands twice a day • Post Infection Support <ul style="list-style-type: none"> ○ Apply topically to infected area at least twice a day

Mechanism of Action

SilverSol Technology[®] is different from all other silvers because of its unique patented technology. Researchers have documented that SilverSol[®] works differently and has multiple modes of action against pathogens.

- SilverSol[®] is metallic silver particle with a thin multivalent silver oxide skin
 - It kills on contact
- SilverSol[®] interferes with viruses
 - Disrupts viral DNA
- SilverSol[®] emits a Resonance
 - Similar to UV light, SilverSol[®] produces a resonate frequency at 890-910 THz, which is antimicrobial
- Wound studies show that SilverSol[®] helps activate stem cells for improved and faster healing

Mechanism of Action

- Silver Oxide
 - Kills on Contact
- Viral Interference
 - Magnetic viral DNA disruption
- Resonance
 - UV light is germicidal at a wavelength of 200-300 nm and a resonant frequency of 890-910 THz
- Stem Cell Activation

With SilverSol[®]

Slide from presentation showing the Mechanisms of action

Test performed by University of Utah showing how SilverSol[®] reduces the inflammation